

# **Simulator Checklist for MD 152**

## **Arterial Stick Arm Simulator**

***Note: Make sure all the following items listed below are in the carrying case/container before re-shelving the simulator***

- 1 Arterial Stick Arm with attached artery tube and red squeezable bulb (pulse) attached to the end
- 1 bottle of Blood Concentrate
- 1 bottle of Manikin Lubricant
- Set of 5 small 3ml Syringes
- Set of 4 replaceable artery tubes
- 1 replaceable skin for arterial stick arm
- 1 replaceable red squeezable bulb (pulse)
- 1 plastic bag containing hardware: screws, bolts, and washers, etc.

### **The Following Items Should Be Inside White Envelope**

- 1 ten-page “Directions for Use” manual
- 1 page sheet titled “Arterial Arm Stick Kit”
- 1 “Medical Definitions for MD 152” list
- 1 “Avoid Damage to Simulator” insert