

Simulator Checklist for MD 132

Lumbar Epidural Injection Trainer

Note: Make sure all the following items listed below are in the carrying case/container before re-shelving the simulator

- 1 Epidural Injection Trainer (includes white skin pad attached to light brown plastic base with removable plastic cover with suction cups attached to it)
- 1 small bottle of simulated Cerebro-Spinal fluid
- 1 packaged replacement 10ml syringe
- Set of 2 small packaged sterile needles
- 1 white plastic container with lid attached to bottom of light brown plastic base containing:
 - a. 1 light blue 10ml “loss of resistance” syringe
 - b. 1 skinny needle with light blue handle
- Set of 4 long clear fluid-filled tubes inside clear plastic sleeve

The Following Items Should Be Inside White Envelope

- 1 Four page Instructions manual
- 1 “Medical Definitions for MD 132” list
- 1 “Avoid Damage to Simulator” insert